




## Trilogy for Mechanical Ventilation Quick Set Up Guide

### ASSEMBLY:








- If the Trilogy has never been used before, it will default to the Full Access Screen. If it is not on the Full Access Screen see below steps:
  - From the Home Screen go to Options.
  - Select Device Options.
  - Set menu access level default to “Full”.

### HOW TO CONFIGURE SETTINGS AND MONITOR A PATIENT:

- On the Home Standby window  select “New Patient” to clear logged patient data and prescriptions, preparing the device for new patient use.
- Select “Other” as the prescription name and acknowledge the use of a bacterial/viral filter with the circuit.
- Tap the  icon to set values for the ventilator. Once on the Setting Screen, select any setting in the selection grid to bring up that setting’s user control in the space below the grid.
- To use as an inpatient vent:
  - Select Circuit: Dual Limb
  - Select Circuit Size: Adult (20-22mm)

- If using an HME and not Passover humidification turn Active Humidification “Off”.
- Choose your Mode of Ventilation by selecting “Mode”:

Trilogy Evo equivalent	Description
A/C-VC	<p><b>Assist Control (Volume Control) mode</b> provides volume-controlled mandatory or assist-control breaths. The set inspiratory time applies to all breaths.</p> <p>If you want to replicate CV mode where the ventilator triggers and cycles all breaths then set the trigger type to OFF.</p>
A/C-PC	<p><b>Assist Control (Pressure Control) mode</b> provides pressure-controlled mandatory or assist-control breaths. The set inspiratory time applies to all breaths. <i>Optional: AVAPS.</i></p> <p>If you want to replicate T mode where the ventilator triggers and cycles all breaths then set the trigger type to OFF.</p>
CPAP	In <b>Continuous Positive Airway Pressure mode</b> , all breaths are spontaneous with the CPAP set pressure delivered in both inhalation and exhalation.
SIMV-PC	<b>Synchronized Intermittent Mandatory Ventilation (Pressure Control) mode</b> is a pressure control mode that provides a mixture of mandatory, assist-control and spontaneous breaths with optional pressure support. It guarantees one mandatory breath in each cycle. The breath rate determines the length of the cycle. <i>Optional: Inspiratory Time min/max for the spontaneous breaths.</i>
SIMV-VC	<b>Synchronized Intermittent Mandatory Ventilation (Volume Control) mode</b> is similar to SIMV-PC, but with volume control.

- After choosing the appropriate mode, make sure AVAPS is turned “Off” in the lower right hand corner of the screen.
  - When settings are complete, select “Accept” to save new values. An unsaved  change indicator is visible until you select “Accept”.
  - To set alarms, select the alarm bell  and set the appropriate alarms. When finished, select “Accept” and “Start Ventilation” to begin therapy.
  - \*\*If a FiO<sub>2</sub> sensor is installed on the back of the Trilogy, from the Home Screen, select  icon then Device Options. From there, select FiO<sub>2</sub> Sensor and turn it “On”\*\*.
- To monitor patient during ventilation, select the  icon to get back to the Home Screen. From the Home Screen, select “Views” in the upper left hand corner of the screen and change view to waveform . The inspiratory phase of the waveforms is color-coded. Orange indicates a ventilator-initiated breath, while blue indicates a patient-initiated breath.
  - To change ventilation modes while using the ventilator or in standby, select the prescription name “Other” to open the prescription list. Select “Add Prescription” then select the name and edit as needed.
  - To put the ventilator in **Standby Mode**, press the On/Off button on the front of the ventilator.